

## Determination of the Refractive Index of $Y_2O_3$ on Glass and Pre-evaporated Substrates by Spectroscopic Ellipsometry

Céline Marchand - Application Scientist - Thin Film Division

$Y_2O_3$  films are good materials for use in many applications such as protection of aluminium and silver mirror coatings, intermediate layers in wide band visible AR coatings and for XeCl laser AR and dielectric mirror designs.

The  $Y_2O_3$  films are hard and usually amorphous with high adhesion to glass, germanium, silicon, zinc sulphide and zinc selenide as well as to metals such as aluminium and silver. In some cases, a very thin layer of yttria can serve as an adherence promoter for multilayer coatings on non-oxide substrates. This material is of medium-index and very transparent over the near-UV range (300 nm) to IR (11  $\mu$ m) region.

Many methods have been used to produce yttria films such as reactive thermal deposition, e-beam deposition, ion assisted deposition (IAD) and magnetron sputtering. As the refractive indices are dependent on the process it is important to know with high accuracy the optical constants and the thicknesses of the films, especially when they are used as high refractive index materials and in combination with higher index materials such as  $TiO_2$  and  $Ta_2O_5$ .

This note describes how the UVISEL Spectroscopic Phase Modulated Ellipsometer was used to characterize the influence of substrate to growth of  $Y_2O_3$ , that is the difference of optical properties of the  $Y_2O_3$  film grown on amorphous glass substrate and the one grown on pre-evaporated layers of  $ZrO_2$ .

### Sample Preparation

The films were prepared by electron beam evaporation in a reactive oxygen atmosphere. The starting materials were grains of  $Y_2O_3$ . The glass substrates were 5 mm thick and were rotated during deposition to improve the uniformity of the films. In the first run two bare glass substrates were put into the chamber and a  $ZrO_2$  film deposited. For the second run one bare glass substrate and the sample from the first run with pre-evaporated  $ZrO_2$  were placed into the chamber and a  $Y_2O_3$  layer was deposited.

Three samples were prepared and characterized as follows.

Sample 1:  $ZrO_2$  monolayer (See Application Note, Ref: SE06),

Sample 2:  $Y_2O_3$  monolayer,

Sample 3:  $Y_2O_3$  on  $ZrO_2$  pre-evaporated.




### Results

The work was performed using the Jobin Yvon UVISEL Visible Spectroscopic Ellipsometer. The Ellipsometric measurements were made at an angle of incidence of  $70^\circ$  across the spectral range 300-830 nm. Both the refractive indexes and thicknesses were extracted from the SE data analysis. The optical constants were determined using the Lorentz Oscillator formula below:

$$n^2 = \epsilon = \epsilon_\infty + \frac{(\epsilon_s - \epsilon_\infty)\omega_t^2}{\omega_t^2 - \omega^2 + i\Gamma_0\omega}$$

## Characterization of $Y_2O_3$ on glass




$Y_2O_3$  is a medium-index and low absorption material. Some amount of index inhomogeneity can appear with increasing layer thickness. The effect can be reduced by providing sufficient oxygen backfill during evaporation. As with the  $ZrO_2$  sample (AN SE-06) a great improvement to the  $\chi^2$  value was found with a two layer model, with the top layer being porous. Using the Effective Medium Approximation (EMA) DeltaPsi2 (DP2) software can determine the degree of porosity compared to the layer on the substrate, and in this case it was found to be 25%.



## Characterization of $Y_2O_3$ on $ZrO_2$

As a first step this sample was characterized using the optical constants found from the previous results. Good values for the films were obtained, but it was found that the results could be improved by the use of a slightly more complicated structure.

Several approaches were followed (interface layer between the two films,...), but the final model that significantly improved the fit, was the optimisation of the  $Y_2O_3$  thicknesses, indices and porosity parameters. In the final results the  $Y_2O_3$  variation in refractive index between two samples was about 0.02.



## Conclusion

The deposition conditions have generated inhomogeneous porous layers that have been observed by the UVISSEL Spectroscopic Phase Modulated Ellipsometer. Moreover, the spectroscopic measurements allow determination of the small discrepancies between the optical constants of a  $Y_2O_3$  film grown at bare substrate and the one grown on a  $ZrO_2$  film.

**HORIBA**JOBIN YVON

(All HORIBA Jobin Yvon companies were formerly known as Jobin Yvon)

Find us at [www.jobinyvon.com](http://www.jobinyvon.com) or telephone:

USA: +1-732 494 8660	France: +33 (0)1 64 54 13 00	Japan: +81 (0)3 3861 8231
Germany: +49 (0)89 462317-0	UK: +44 (0) 20 8204 8142	Italy: +39 02 57603050
China: +86 (0)10 6849 2216	Other Countries: +33 (0)1 64 54 13 00	