





PTI QuantaMaster™ Series

Modular Research Fluorometers for Lifetime and Steady State Measurements



Most sensitive and flexible
fluorometer

Unique PTI QuantaMaster Benefits

- *Highest guaranteed sensitivity specification*
- *World class software for all steady state and lifetimes needs*
- *Extended wavelength range with triple grating monochromators*
- *Multiple light sources and detectors with dual entrance & dual exit monochromators*
- *Excellent stray light rejection with large single or double additive, coma corrected, monochromators*
- *World class TCSPC lifetime enhancements, up to 100 MHz*
- *NIR steady state and phosphorescence lifetime detection to 5,500 nm*





“The new PTI QuantaMaster™ 8000 series modular research fluorometers from HORIBA Scientific offer the world’s highest guaranteed sensitivity specification, plus many unique benefits.”


PTI QuantaMaster 8000 Series

The PTI QuantaMaster™ series of modular research grade spectrofluorometers are multidimensional systems for photoluminescence measurements. The foundation of a fluorescence spectroscopy laboratory is built on steady state intensity measurements such as wavelength scans, time-based experiments, synchronous scans and polarization. The PTI QuantaMaster series ensures you get the best possible results for all these measurements with high sensitivity, spectral resolution and stray light rejection. This level of sensitivity is achieved using a unique xenon illuminator, providing safety, cost, and energy consumption benefits not found amongst competitor companies. These conditions make the PTI QuantaMaster system ideal for the highest demands of fluorescence research.


The PTI QuantaMaster system is adaptable to every research need, with additions such as TCSPC fluorescence lifetimes, upconversion lasers and phosphorescence detection up to 5,500 nm. Using an optional pulsed light source allows for not only spectral and kinetic fluorescence and phosphorescence measurements, but also the measurement of lifetimes in the picosecond to seconds range. This addition is especially beneficial when using fluorescent probes prone to photobleaching, and when characterizing inorganic material with longer lifetimes. The modular design of the PTI QuantaMaster series ensures that your system can be easily adapted to your growing research needs.

Ultimate Sensitivity

The industry standard for sensitivity of a fluorometer is the signal to noise ratio calculated from a water Raman spectrum. Using this standardized test (see our technical note "SN Determination for PTI QuantaMaster Fluorometer,") our signal to noise ratio is the highest sensitivity in the world. The extreme sensitivity of the PTI QuantaMaster fluorometer is achieved with the lowest wattage lamp in the industry. This is a result of the intelligent engineering of the unique PTI PowerArc™ arc lamp illuminator featuring an ellipsoidal reflector with the highest possible light gathering efficiency of 67%, and focusing the light in a tight spot at the monochromator slit. As a result, the standard 75W Xe delivers light to the sample more efficiently than higher power lamps featured by other instruments. This reduces energy waste and excessive heat generation by an overpowered light source, not to mention cost, while exceeding the sensitivity of all competitors' designs. Another contribution to the high sensitivity of the PTI QuantaMaster comes from the asymmetrical, aberration-corrected monochromator, optimized for the best light throughput and stray light rejection.




Water Raman spectrum of the PTI QuantaMaster 8000, resulting in a signal to noise ratio of >30,000:1! Experimental conditions: $\lambda_{\text{exc}} = 350 \text{ nm}$, $\Delta\lambda_{\text{exc}} = \Delta\lambda_{\text{em}} = 5 \text{ nm}$, $\text{int} = 1 \text{ s}$.



Superior optical design of the PTI QuantaMaster 75W Xe lamp illuminator is a key factor in PTI QuantaMaster 8000 sensitivity

Want a higher wattage light source?

No problem. We offer a 450W lamp that retains the same system sensitivity.



Ultimate Stray Light Rejection


With Single or Optional Double Additive Monochromators

Suppression of stray light is one of the most critical factors when measuring highly scattering, or low quantum yield samples. Every PTI QuantaMaster 8000 series spectrofluorometer is custom made with the highest quality optics to insure the lowest amount of scatter. The standard 300 mm focal length asymmetrical Czerny-Turner monochromators are coma-corrected and individually optimized for purpose as either excitation or emission monochromators, ensuring the lowest amount of stray light contamination for the best detection of the true fluorescence signal. These monochromators boast an impressively high stray light rejection of 1×10^{-5} in a single excitation monochromator configuration. For more sensitivity and higher performance, the PTI QuantaMaster 8000 can also be configured with double additive 300 mm focal length monochromators, improving the stray light rejection to 1×10^{-10} . The PTI QuantaMaster 8000 also offers an optional order sorting filter wheel for rejection of second order signals from large spectral scans. This ultimate stray light performance was motivated by an increasing demand for photoluminescence spectrometers in materials science, where strongly scattering samples, such as powders, wafers and films, are routinely used. Very low stray light performance will also benefit researchers working in biological, biomedical and environmental areas where cell suspensions, protein and biomembrane solutions, or soil samples, generate high levels of scattered light.


Resolution

Resolution is of utmost importance to photoluminescence research. High quality resolution can reveal detailed spectral features which is indispensable for applications in materials science and analytical chemistry. Resolution is the key to detecting very narrow lines, which is necessary to study fine interactions in inorganic materials and crystals. The PTI QuantaMaster 8000 yields high quality resolution due to innovative optical design, and very minimal optical aberrations.

The PTI QuantaMaster 8000 spectrofluorometers use a precision-driven asymmetrical, 300 mm focal length Czerny-Turner monochromator with a motorized triple grating turret and motorized flipping mirrors. More than 30 different gratings are available. Due to the combination of the computer-controlled motor with micro-stepping resolution and available grating selection, it is possible to achieve a minimum 0.022 nm step size. This means that within the UV and Vis spectral regions, you can resolve spectral features well below 0.1 nm.



Raman spectrum of CCl₄ using PTI QuantaMaster 8000 with single excitation & single emission monochromators. Well-resolved peaks and no contamination from the Rayleigh scattering demonstrates the excellent stray light suppression. Experimental conditions: $\lambda_{\text{ex}} = 349 \text{ nm}$, $\Delta\lambda_{\text{ex}} = 0.7 \text{ nm}$, $\Delta\lambda_{\text{em}} = 0.7 \text{ nm}$, step size=0.05 nm, integration time = 1 s



Emission scan of dysprosium-doped YAG crystal measured at 78K in software controlled LN cryostat illustrating excellent resolution of narrow spectral lines attained at low temperature. Experimental conditions: $\lambda_{\text{ex}} = 353 \text{ nm}$, $\Delta\lambda_{\text{ex}} = 5 \text{ nm}$, $\Delta\lambda_{\text{em}} = 0.1 \text{ nm}$, step size=0.022 nm, integration time = 1 s

Spectral Range and Signal Detection


For most applications, the typical detector employed is a PhotoMultiplier Tube (PMT). The standard PTI QuantaMaster 8000 configuration features a highly sensitive PMT, with the option of photon counting, analog, TCSPC and SSTD detector modes. The PTI QuantaMaster Series offers you the ability to customize the system to meet your applications needs. Digital detection, or photon counting, offers the highest sensitivity as it records single photon events. The analog detection measures the current that is generated on the PMT anode, and provides for additional detection gain control. This greatly enhances the dynamic range of the instrument, especially for higher intensity signals. For NIR and IR applications, we also offer specialized PMTs and solid state detectors, such as InGaAs, PbS and InSb diode detectors that are capable of detecting out to 5500 nm. Most of these detectors can be used with pulsed light sources for time-resolved photoluminescence. Multiple detectors can be used with a single instrument: a single monochromator will accept two, and a double monochromator, up to three detectors. The selection is done by computer-controlled steering mirrors which direct the emitted light to a selected detector. A triple, motorized grating turret ensures good light efficiency for any detector range.



Excitation and Emission Correction

All light sources emit light that is not of equal intensity across the output spectrum, and this can lead to errors in the measurement of an excitation spectrum. The raw data must then be corrected for this discrepancy. The PTI QuantaMaster 8000 utilizes a reference diode detector that has been calibrated and installed at the factory. Excitation correction is performed in real time. During an experiment, part of the excitation beam is diverted prior to reaching the sample. This fraction of photons is measured and then the reference detector provides a correction that is independent of the excitation source characteristics, or any temporal fluctuation of the lamp intensity, thus ensuring excellent stability of the signal.

A similar phenomenon exists for emission spectra. Since the detection efficiency of the optics, gratings, mirrors and detector is not equivalent at all wavelengths, some type of correction must be performed to account for these variations. Typically, the emission channel is calibrated at the factory with a known light source, such as a NIST-traceable standard. This information is used to construct a correction file, which is then stored locally on your computer. Multiplication of the raw data by this correction file yields the true corrected emission spectrum. This correction can be performed in real-time, or can be recalled in later analysis of the raw data and applied in the easy-to-use PTI FelixGX software.




Raw and corrected pyrene excitation and emission spectra with excimer peak present around 475 nm. Corrected data shown in red.

Specifications

The following specifications are for the standard PTI QuantaMaster 8075-11 and 8450-11 systems. Options and upgrades may be available upon request.

Signal to Noise Ratio	>30,000:1 RMS, 15,000:1 FSD
Data Acquisition Rate	1,000,000 points/sec. to 1 point/1000 sec
Inputs	4 photon counting (TTL); 4 analog (+/- 10 volts); 1 analog reference channel (+/- 10 volts); 2 TTL
Outputs	2 analog (+/- 10 volts); 2 TTL
Emission Range	185 nm to 900 nm (optional to 5,500 nm)
Light Source	High efficiency "ECO" friendly continuous 75 W Xenon arc lamp (Optional 450 W Xenon)
Monochromator	300 mm, triple grating, coma-aberration corrected, asymmetrical, excitation or emission optimized, Czerny-Turner design
Slits	Computer controlled, continuously adjustable
Excitation Grating	1200 line/mm 300 nm blaze, (Up to two optional gratings can be ordered)
Emission Grating	1200 line/mm 400 nm blaze, (Up to two optional gratings can be ordered)
Wavelength Accuracy	+/- 0.3 nm
Minimum Step Size	0.022 nm
Standard Detection	Multimode: Photon Counting, 3 analog (fast, medium, slow response), direct and Single-Shot Transient Digitizer (SSTD) mode, and Time-Correlated Single Photon Counting (TCSPC)
System Control	Computer interface with FelixGX spectroscopy software
Lifetime Range	5 ps to seconds with appropriate time-resolved accessories
Dimensions (W x D x H)	36 in x 28 in x 12 in




TCSPC Lifetime Measurements


The PTI QuantaMaster 8000 series can be easily enhanced with TCSPC fluorescence lifetime capabilities. Utilizing world class TCSPC sources, electronics and detectors developed by our IBH group, the PTI QuantaMaster provides the ultimate in speed, versatility and performance. The standard QuantaMaster PMT can be used for these additional TCSPC measurements, or you can add a dedicated TCSPC detector with enhanced performance or extended NIR detection. All steady state and time-resolved control, acquisition and analyses are handled by FelixGX software.

HORIBA TCSPC Benefits

- 40 years of experience in TCSPC innovation
- Industry-leading true 100 MHz system operation allows for millisecond acquisition times
- TCSPC lifetime measurements from 5 ps to seconds
- Full control over TCSPC and steady state acquisitions with single FelixGX software package
- Measure TCSPC lifetimes, time-resolved anisotropy and TRES (Time-Resolved Emission Spectra)
- Select from our catalog of over 60 state-of-the-art compact pulsed LEDs and Laser Diodes for virtually any application
- For unsurpassed versatility, choose a picosecond supercontinuum laser with HORIBA's proprietary Frequency Doubler – a powerful tool for time-resolved protein studies
- PowerFit-10 decay analysis package with multiple fitting models including a unique Maximum Entropy Method (MEM) lifetime distribution program



Fluorescence decay of CdSe QDots measured with TCSPC option of QuantaMaster. The MEM lifetime distribution (inset) reveals size heterogeneity of QDots




PTI QuantaMaster 8000 with DeltaTime TCSPC option and frequency doubled Supercontinuum laser is a perfect choice for intrinsic protein TR fluorescence

Two Systems Can Be Better Than One!


We also offer very affordable stand-alone TCSPC systems. You can increase your lab's throughput by having a dedicated steady state fluorometer and a dedicated TCSPC system operating at the same time, for almost the same price as adding TCSPC to the PTI QuantaMaster.

Near-Infrared Spectrofluorometry

Near-infrared (NIR) spectrofluorometry has emerged as a valuable analytical technique, especially in the fields of material research, nanotechnology, chemistry, and photomedicine. Powerful and diverse NIR capabilities are available from HORIBA as either a stand-alone research grade fluorometer, or as an upgrade to our UV-Vis steady state spectrofluorometers. There are different configurations to adapt to any research needs.



PL excitation and emission spectra of Er^{3+} doped glass measured with the PTI QuantaMaster 8000 equipped with the NIR TE-cooled InGaAs detector.



Er decays: PL decays of Er^{3+} doped glass at different emission peaks measured with the PTI QuantaMaster 8000 system equipped with the NIR TE-cooled InGaAs detector and pulsed nitrogen/dye laser excitation.

NIR-PMT Based Detectors

These detectors offer the ultimate in sensitivity and can be used for steady state and TCSPC measurements.

Available with four NIR PMTs for maximum spectral range coverage:

- 300 nm–1,400 nm, LN-cooled
- 950 nm–1,400 nm, TE-cooled
- 300 nm–1,700 nm, LN-cooled
- 950 nm–1,700 nm, TE-cooled

Solid-State Photodiode-Based NIR Detectors

There are a variety of photo diodes available, with TE- or LN-cooling and a lock-in amplifier and chopper for enhanced sensitivity.

- InGaAs: 500 nm–1,700 nm (up to 2400 nm detector available upon request)
- PbS: 1,000 nm–3,200 nm
- InSb: 1,500 nm–5,500 nm

Additions for NIR Lifetime Measurements to 5,500 nm!

All of the detectors listed above can be used in a Single Shot Transient Digitizer (SSTD) mode for luminescence lifetime measurement capability in NIR to measure lifetimes from 1 μs to hundreds of ms. SSTD is extremely fast and offers outstanding signal to noise, using:

- Variable high rep rate pulsed xenon lamp option for phosphorescence lifetime (NIR-TR-10)
- Pulsed nitrogen and dye laser (NIR-TR-20)
- 3rd party pulsed Q-switched lasers
- TCSPC lifetime add-on with supercontinuum laser, laser diodes and LEDs (for NIR-PMT based systems)


Flexibility

Enhanced Measurements



Unique NIR Solutions

A PTI QuantaMaster 8000 can be equipped with multiple illuminators and detectors to cover the widest spectral range for both steady state spectra and fluorescence and phosphorescence lifetimes. Consider the following configuration:

- Double emission monochromator with R928 PMT, InGaAs & InSb detectors cover 250 to 5,500 nm
- Continuous xenon lamp for steady state spectra
- 20 Hz Q-switched/OPO Opolette laser for tunable excitation from 210 to 2,200 nm
- Frequency doubled nitrogen pumped dye laser for tunable excitation from 250 to 990 nm
- Steady state spectra from 250 to 5,500 nm
- Single Shot Transient Digitizer (SSTD) for phosphorescence decays over entire range from 250 to 5.500 nm




Ho ion emission and excitation spectra and PL lifetimes in NIR measured with cw Xe lamp and pulsed N₂/dye laser excitation with the same InGaAs detector






Unique NIR Solutions


Flexibility



Spectra and PL lifetimes in NIR of Er-doped glass measured with continuous Xe lamp and pulsed N2/dye laser excitation with the same InGaAs detector



Nd decay with PbS detector - deconvolution required




Modularity

Many Configurations to Fit your Specific Needs, and Modularity to Grow

The PTI QuantaMaster 8000 series features an open architecture design that provides the ultimate in versatility to adapt to any future fluorescence application needs. You can optimize the initial configuration by choosing the light source, gratings, and PMTs, as well as a wide array of available accessories. The number of available configurations is virtually limitless!

The PTI QuantaMaster 8000's universal QuadraCentric™ sample compartment has a spacious design that provides accessibility, and can accommodate a wide selection of sample accessories. Choose from sample temperature controllers to various holders for solids, liquids and powders, dewars, integrating spheres and many other options. See the Accessories page for more details.




PTI QuantaMaster 8075-11 with a second emission channel

QM-8075-22 fluorometer with double additive excitation and double additive emission monochromators.

PTI QuantaMaster 8075-11 with TCSPC lifetimes with DeltaDiodes or NanoLEDs and detection

The open architecture design also allows for application and methodology changes. As your application needs grow, so can your PTI QuantaMaster. For example, if you develop a need to measure dynamic anisotropy, you can add a second emission channel and a set of polarizers. If you want to complement your steady state data with lifetime measurements, you can do so by adding a laser or LED-based excitation to your initial configuration. If you are interested in intracellular Ca^{2+} after completing initial Fura-2 studies, you may decide you would like to start imaging the events. The system can be easily coupled with any fluorescence microscope. Whether you choose to add NIR detection or a second excitation source, the possible configurations are endless.




PTI QuantaMaster 8075-11 upgraded to fluorescence microscopy with an additional PMT detector equipped with an eyepiece aperture

PTI QuantaMaster 8075-11 with a pulsed light source and a gated detector for phosphorescence or lanthanide emission

PTI QuantaMaster 8075-11 coupled to a microscope with detection into the existing emission monochromator

FelixGX Software

PTI QuantaMaster 8000 fluorometers come with our integrated FelixGX software to control both the instrument and accessories. Designed to be complete and easy to use, FelixGX includes analytical functions for trace manipulation, and spectral and kinetic analysis. Connected with a USB interface, FelixGX provides a full set of data acquisition protocols, and controls the hardware for all system configurations and operating modes.



Software

FelixGX Controls

Hardware Controls

- Single or double monochromators
- Triple grating turret
- Flipping mirrors switching between different light sources and detectors
- Motorized slits
- Motorized polarizers
- Motorized multiple sample holders
- Excitation correction detector (Xcorr)
- Temperature control Peltier devices
- Cryostat
- Gain control of PMT detector
- Switching from digital to analog mode
- External devices such as stopped flow and titrator
- Pulsed light sources
- Scanning of wavelength-tunable OPO lasers
- TCSPC electronics
- Electroluminescence and photovoltaic accessories
- Supercontinuum lasers


Acquisition Modes

PTI FelixGX provides several acquisition modes for spectral and kinetic measurements:

- Excitation and Emission spectral scans with user control of integration time, monochromator step, speed and wavelengths
- TCSPC
- Time-based scan with user defined macro-time duration and integration time
- Spectral and time-based polarization scans with full control of motorized polarizers and automated measurement of G-factor, and sample background for all polarizer orientations
- Simultaneous multi-dye measurements with pre-defined library of common fluorescence dyes or customer-defined dyes
- Synchronous excitation/emission scan
- Excitation and emission ratio fluorescence
- Phosphorescence decay and time-resolved excitation and emission spectra using Single-Shot Transient Digitizer (SSTD)

Macro Capabilities

FelixGX comes equipped with Macro capabilities to allow for automated measurements. Choose from a list of actions to make a chain of commands, or set up a loop function to eliminate the need to constantly change the acquisition settings. Set up the automation job and simply walk away, letting FelixGX execute your instructions.




FelixGX Analytical Capabilities

Trace Manipulation

FelixGX also provides an extensive set of math functions that can be used for processing and manipulation of acquired data traces:

- Antilog
- Average
- Distribution average
- XY Combine
- Differentiate
- Integrate
- Linear Fit
- Linear Scale
- Logarithm
- Normalize
- Reciprocal
- Smooth
- Truncate
- Baseline
- Merge Traces
- Peak Finder




Demonstration of Peak Finder function with an emission spectrum

Kinetic Data

Fluorescence and phosphorescence decays can be analyzed with the TCSPC lifetime analysis package which includes:

- One to four exponentials
- Multi one to four exponentials
- Global one to four exponentials
- Anisotropy decays
- Exponential Series Method (ESM) lifetime distribution
- Maximum Entropy Method (MEM) lifetime distribution
- Micelle kinetics
- Stretched exponential
- DAS (Decay Analysis Software)/TRES



TCSPC output graph with residuals and autocorrelation.

Advanced Calculators

FelixGX also offers a special set of software functions, such as quantum yield, absorption, FRET and color coordinates calculators, as well as the software that calculates structural parameters for single-walled carbon nanotubes. These are very convenient additions to some accessories, such as the integrating sphere or absorption accessory, and are also indispensable for some fluorescence applications, such as intermolecular interactions (FRET) and materials characterization.

Absorption Calculator

Absorption measurements are complementary to fluorescence. They are necessary for fluorescence quantum yield determination and are an easy and convenient way to check the fluorophore concentration. You can compare the absorption and excitation spectra to draw conclusions about the purity of the sample. Using the built-in absorption calculator with an absorption accessory will greatly enhance the capabilities of your PTI QuantaMaster 8000 fluorometer.

The screenshot shows the 'Absorption Accessory Calculator' window. It has two main modes: 'Spectral' (selected) and 'Timebased'. In 'Spectral Mode', users can input 'Low X' (86.95) and 'High X' (194.05). They can also select 'I₀(λ)' and 'I(λ)' from a dropdown menu, both currently set to 'D1 450:500-900 ex. open em. 1nm'. There are buttons for 'Calculate Absorbance' and 'Calculate Transmittance'. The 'Timebased Mode' section includes input fields for 'I₀(t)', '<I(t)>', 'I(t)', and '<I(t)>', each with a 'Calculate' button. At the bottom, there are fields for 'Absorbance' and 'Transmittance'.

Quantum Yield

Quantum yield is one of the most important parameters that characterize photoluminescence of materials. FelixGX incorporates a quantum yield calculator which, when coupled with an integrating sphere, allows you to calculate the quantum yield with ease.

The screenshot shows the 'Quantum Yield' calculator window. It is divided into 'Emission Traces' and 'Excitation Traces' sections. 'Emission Traces' includes 'I_{em}(λ)' (A1 325:350-900 SPL) and 'I_{ref 2}(λ)' (none), with an 'Execute 1' button. 'Excitation Traces' includes 'I_{ex}(λ)' (A1 325:280-400(2) b) and 'I_{ref 1}(λ)' (A1 325:280-400 SPL), with an 'Execute 2' button. On the right, 'Integrals' shows values: 5192.738, 0, and Diff: 5192.738. Below, 'Range' is set with 'Low X: 280.00000' and 'High X: 399.45370'. 'Scaling' shows a 'Scaling Factor' of 0.0002. At the bottom, 'Trace Pair Adjustment' has 'Excitation' selected. The 'Quantum Yield Result' section shows 'Quantum Yield: 0.000898629' and a 'Calculate' button.

FRET



The screenshot shows the 'FRET - Determine R0' window. At the top, it displays the Förster equation:
$$R_0 = 0.2108 \sqrt{\kappa^2 \Phi_D n^2 \frac{\epsilon(\lambda_{max})}{E_A(\lambda_{max})} \frac{\int I_D(\lambda) E_A(\lambda) \lambda^4 d\lambda}{\int I_D(\lambda) d\lambda}}$$
 Below the equation, 'Data Curves' are set to 'Donor Emission' (D1 200-460:465) and 'Acceptor Absorption' (D1 200-460:465 [COR]). The wavelength λ_{max} is set to 435. 'Parameters' include $\kappa^2 = 0.6666666$, $n = 1.3333333$, $\Phi_D = 1$, and $\epsilon(\lambda_{max}) = 20000$. A 'Set To Default' button is present. At the bottom, 'Förster distance (Å)' is calculated as 42.319, with a 'Calculate R0' button and a 'Close' button.

The FRET (Förster Resonance Energy Transfer) technique provides information about molecular distances, interactions in macromolecular systems, binding, diffusion, sensing, etc. FRET happens when an excited donor molecule transfers its energy to an acceptor in the ground state. FRET is essentially a molecular ruler, where distances are scaled with the Förster critical radius R_0 , which is a unique parameter for a given donor-acceptor (D-A) pair, defined by spectroscopic parameters of the pair and their environment. Once the R_0 is known and the FRET efficiency is determined experimentally, the D-A distance and the FRET rate constant can be calculated. From spectra or TCSPC data, FelixGX provides an easy and convenient way of calculating all relevant FRET parameters, including R_0 .

Advanced Calculators



Single-Walled Carbon Nanotube Calculator

Carbon nanotubes can be characterized using the specially-designed NanoCal within FelixGX. NanoCal analyzes 3-D ExEm spectral maps and returns structural parameters such as the nanotube radius and the chiral angle. Combining this easy-to-use software with PTI QuantaMaster 8000 NIR options allows for full characterization of SWCNTs.



Color Coordinate Calculator


In many applications, such as phosphors for screen displays, multi-color LEDs, fluorescent additives to consumer products, etc., there is a need to quantify a visual perception of color. FelixGX provides a Color Coordinate Calculator based on two widely accepted standards introduced by the International Commission on Illumination, CIE 1931 and CIE 1976. The CIE 1931 uses x,y chromaticity coordinates where each x, y pair corresponds to a unique color within the colored shape. The CIE 1976 uses a system with more uniform perceptual chromaticity to define the color space using u, v coordinates. Upon highlighting a spectral trace and clicking on CIE 1931 and CIE 1976 Color Coordinates, FelixGX will display both CIE pairs.



Applications and Examples

Steady State Anisotropy


Both photon absorption and photon emission are correlated in space with the transition dipole moment vector of the molecule. Therefore, a measurement of the fluorescence polarization of the emitted light can yield information about the rotational mobility of the molecule under investigation. The rotational mobility of a macromolecule, such as protein or DNA, depends on its size, conformation and viscosity of the medium. Fluorescence anisotropy measurements provide an easy and powerful tool to study conformational transitions, such as protein folding and unfolding induced by temperature, pH changes, and drug or ligand binding. For fast and convenient anisotropy measurements, dual emission configurations are available to allow simultaneous determinations of vertically and horizontally polarized fluorescence signals. A software controlled rapid temperature change Peltier unit is a valuable option for anisotropy measurements.




Temperature-induced unfolding of bovine serum albumin (BSA) in PBS (pH=7.4) monitored by fluorescence anisotropy with dual emission channels using a rapid temperature change Peltier option.

Total Luminescence Spectroscopy (TLS)

The powerful FelixGX software, with its user-friendly acquisition macro programming capability, and the rapid scanning performance of the PTI QuantaMaster 8000, make it easy to create automated acquisition protocols for measuring emission spectra at varying excitation wavelengths, and creating a 3D EEM characterization of a fluorescing sample. Such measurements enable the user to fully characterize spectrally complex samples rapidly and automatically, saving valuable time. The TLS technique is used in various analytical applications of photoluminescence spectroscopy. It is especially useful for detecting and identifying Polycyclic Aromatic Hydrocarbons (PAHs) in environmental samples, as well as in food science to test for contaminants, or assess foodstuff deterioration. 3D mapping can then be used to demonstrate the naturally occurring fluorescent components.



Rapid automated Ex/Em Matrix scan of p-terphenyl/anthracene mixture in a 3-D representation.




TLS is an efficient method to assess the quality of beer. Beers exhibit complex intrinsic fluorescence with contributions from tyrosine, tryptophan (300-400 nm) and flavins (400-500 nm). The diminished contribution from flavin fluorescence is indicative of deteriorating taste and quality of beer (e.g. 'skunked' beer).

Applications and Examples

Synchronous Fluorescence Spectroscopy (SFS)


Synchronous Fluorescence Spectroscopy involves scanning the excitation and emission monochromators simultaneously at identical scan rates, with a fixed offset between the two wavelength ranges. It offers much higher spectral selectivity than the conventional emission and excitation scans, reduces light scattering, and improves resolution. SFS is a powerful analytical technique that enables simultaneous determination of multiple components in the mixture. It has been used in detecting carcinogenic Polycyclic Aromatic Hydrocarbons (PAHs) in food and environmental samples.




The data represents a mixture of three organic hydrocarbons: p-terphenyl, anthracene, and perylene. The ordinary emission scan does not reveal the complexity or identities of the mixture. On the other hand, the synchronous scan clearly shows 3 narrow emission peaks located at the emission maxima of the respective compounds, making it possible to identify the mixture components.

Ratiometric Measurements for Intracellular Ions

Excitation-shifted probes such as Fura-2 and BCECF are often used in determining intracellular calcium concentration and pH, respectively. These probes exhibit an excitation shift upon binding calcium (Fura-2) or protonation (BCECF). In these experiments, the excitation monochromator automatically alternates between two excitation wavelengths corresponding to the free and ion-bound probe. The ratio of the two signals is then measured. Pre-configured look-up tables transform the measured ratio into ion concentrations or pH. Similar measurements can be done for emission shifted probes such as Indo-1 and carboxy-SNARF.



Fura-2 titration with Ca²⁺ ions monitored via excitation spectra.




Freshly isolated rabbit atrium was stimulated with 90 mM KCL. 2 μ M nicardipine was added and returned to normal Tyrode medium, followed by the addition of 10 mM caffeine. The Fura-2 excitation ratio signal follows the kinetics of free Ca²⁺ in the tissue.

Time Based Measurements

Probably one of the most common experiments, time based measurements, are useful for many applications such as enzymatic activity assays, ion activity in cells, titration studies, protein-protein and protein-drug interactions, anisotropy measurements, and chemical kinetics. The measurements involve monitoring the fluorescence intensity at fixed excitation, and emission wavelengths as a function of time. The PTI QuantaMaster 8000 series can do kinetic measurements on a time scale ranging from microseconds to hours or days. The use of the excitation correction unit (Xcorr) greatly improves the signal stability by eliminating any light source intensity fluctuations and drift over time.


For best results, time-based kinetic experiments should be conducted at a controlled temperature. Therefore our Peltier-based rapid temperature controlled cuvette holders, K-155-C or K-157-C, are recommended. If very fast reaction kinetics are studied, a stopped-flow accessory, K-161-B, will be a useful addition.



BSA protein unfolding induced by detergent SDS monitored in time-based mode.

Förster Resonance Energy Transfer (FRET)

FRET is a popular technique used to study binding, conformational changes, dissociation and other types of molecular interactions. Applications of FRET are especially common in biomedical research involving protein-protein, protein-nucleic acid interactions, protein folding/unfolding, nucleic acid hybridization, membrane fusion and many others. There is also a variety of immunoassays based on FRET. The FRET phenomenon occurs between an excited donor (D) molecule and a ground-state acceptor (A) molecule over a range of distances, typically 10-100 Å. It is a nonradiative process, meaning no photon is emitted or absorbed during the energy exchange. The efficiency of FRET is strongly dependent on the D-A distance and is characterized by the Förster critical radius R_0 , a unique parameter for each D-A pair. Once R_0 is known, the D-A pair can be used as a molecular ruler to determine the distance, or monitor distance changes between sites labeled by D and A. Since FRET is mostly used to study biological systems, where concentrations are often low, and samples can be highly scattering, the PTI QuantaMaster 8000 is an ideal fluorometer for this application due to its high sensitivity and excellent stray light rejection. It is also easy to upgrade to a lifetime option, which can be very beneficial for verification of the FRET mechanism. The PTI QuantaMaster 8000 series will also help you take advantage of this technology with the built-in FelixGX FRET Calculator.



Titration monitored by FRET between Alexa-BSA complex and a Bodipy-labeled fatty acid.

Applications and Examples

Automated Temperature Control


Sample temperature plays a critical role in all types of luminescence measurements. For example, when the emission anisotropy is measured, the viscosity will change as a function of the temperature affecting the rotational motion of the fluorophore. The temperature control can be critical for fluorescence quantum yield determination, or any quantitative intensity measurements since the nonradiative deactivation is strongly temperature dependent. Temperature control is essential in fluorescence studies of proteins as it affects thermal stability of proteins, and their folding and unfolding characteristics.

Solid samples, such as doped crystals, glasses, ceramics, and organic molecules deposited on surfaces will exhibit narrowing of spectral lines when cooled to low cryogenic temperatures, thus allowing study of fine interactions. Organic molecules will usually exhibit phosphorescence when cooled to sufficiently low temperatures.

The PTI QuantaMaster 8000 series comes standard with a thermostatable cuvette holder where the plumbing is already in place for temperature control utilizing a circulating water bath. If your research requires more precise or extreme temperature control, additional sample holders are available, including software controlled Peltier-based variable temperature cuvette holders (single or 4 position) and a liquid nitrogen cryostat. Programmable spectral scans at automatically varying temperatures and temperature ramping experiments are available.

Temperature control is critical in applications, such as:


- Temperature dependent quantum yields
- Quantitative intensity measurements
- Activation energies of photophysical processes
- Protein folding and unfolding
- Nucleic acid melting profiles
- Thermodynamic parameters of binding reactions
- Membrane fluidity and permeability studies
- Fluorescence measurements of live cells
- Enzyme kinetics




Fully automated temperature mapping emission scans of coronene deposited on silica measured in liquid nitrogen cryostat. As the temperature is lowered, the phosphorescence spectrum begins to appear and intensity increases.

Phosphorescence with a Pulsed Light Source

A pulsed light source and the ability to integrate the signal at user-selectable time delays are indispensable tools in discriminating spectra based on the lifetime of the respective excited state. Fluorescence emission happens on the picosecond to nanosecond time scale, while phosphorescence occurs on the microsecond to second time scale. By varying the temporal position and the width of the signal detection gate, one can selectively detect fluorescence and phosphorescence spectra as attested by phenanthrene spectra on the accompanying figure. Here, the emission of phenanthrene in a frozen glass was measured with gradually increased time delay of the detection gate to diminish contribution of fluorescence. However, the true potential of this technique can be seen in the case of Room Temperature Phosphorescence (RTP) of RNase T1 tryptophan, where the signal was extracted by gating out the overwhelming Trp fluorescence—a task impossible with a continuous excitation source. Conveniently, the same instrument can be used to measure phosphorescence decay of this extremely weak emission by using the Single-Shot Transient Digitizer (SSTD) function of the PTI QuantaMaster interface.




Phenanthrene at 77 K utilizing a cold finger nitrogen dewar accessory. Fluorescence and phosphorescence spectra measured while increasing the delay time (at 2 μ s increments) for signal integration.




Discrimination between strong fluorescence and weak Room Temperature Phosphorescence (RTP) from RNase T1 tryptophan by varying the temporal position and widths of the signal detection gate on a PTI QuantaMaster 8000 equipped with a pulsed Xe lamp and gated detector for signal integration.

Bio and Chemiluminescence




The unsurpassed sensitivity of the PTI QuantaMaster 8000 detection makes it a very capable instrument for measuring extremely weak chemiluminescence emission, as illustrated by the cytochrome C/hydrogen peroxide experiment.




Phosphorescence decay of a weakly emitting RNase T1 tryptophan signal using the same instrument.

Semiconductors Research



Due to its dedicated accessories, such as a well-designed solid sample holder and excellent stray light rejection characteristics, the PTI QuantaMaster 8000 is an excellent choice for semiconductors research. Here, clean spectra from strongly scattering ZnO samples were measured with the PTI QuantaMaster 8000 equipped with a double excitation monochromator.

Fluorescence Upconversion




Upconversion phenomena in lanthanide-doped glasses and powders has been extensively studied in recent years. They are of interest due to a demand for compact and efficient lasers and amplifiers for optical communications, especially operating in blue-green and orange spectral ranges. The data here illustrates the PL spectra of upconverting nanoparticles NaYF₄: Yb,Er in water solution, measured with the UPCONV-980 upconversion kit for the PTI QuantaMaster 8000, consisting of a 980 nm CW diode laser and an integrating sphere. The upconversion setup allows for simple determination of the luminescence upconversion quantum yields due to an efficient and ergonomically-designed integrating sphere (shown with the top removed) with easy access to the sample area.




Optional compact integrating sphere fits inside the PTI QuantaMaster 8000 sample compartment.

Electroluminescence and Photovoltaic Measurements





The flexibility of the modular design makes it easy to utilize the PTI QuantaMaster 8000 for more specialized applications, such as electroluminescence or photovoltaic measurement. Here, the figure shows an electrical response of a photovoltaic cell illuminated with the PTI QuantaMaster excitation monochromator equipped with an NIR grating. The electrical signal from the cell is fed directly to one of the analog inputs of our versatile interface, and the powerful FelixGX software takes care of rest!



Emission generated by applying different voltages (6V and 9V) to a thin film electroluminescent sample using the PTI QuantaMaster 8000 with the electroluminescence accessory.

Quantum Yield




Quantum Yield determination of Nd³⁺ doped glass in NIR with the integrating sphere and InGaAs detector. The measurement requires high signal stability and precise emission corrections. The QY experiment involves emission scanning over the excitation peak, which is usually significantly higher than the emission spectrum. Since the absorbance of the sample is very low, an excellent signal stability, high dynamic range and a linear behavior of the detector are of utmost importance for accurate QY determination. The graph on the right shows the expanded excitation peak with, and without, the sample. Capturing the difference of the two signals is the key to accuracy. The triplicate experiment showed excellent reproducibility resulting in QY = 0.567 ± 0.017.

Applications and Examples: NIR Fluorescence Spectroscopy

The applications and interest in NIR photoluminescence have been growing rapidly in recent years. This trend is spearheaded mostly by extensive research in nanotechnology and materials science. NIR-emitting nanoparticles, lanthanide doped glasses and ceramics used in developing new laser media and photonic devices, single-walled carbon nanotubes, semiconductor and electroluminescent systems are only a few dominant applications. There is also a considerable research effort in the optical fiber communication industry to develop infrared molecular amplifiers for the transmittance window at 1550 nm. In biomedical areas, there is a trend of using NIR-emitting nanoparticles as luminescent markers due to the fact that the light scattering, a notorious problem in UV-Vis fluorescence measurements, is greatly reduced as the wavelength increases. Less interference means better signal to noise with strongly scattering biological samples. NIR light can penetrate tissue at a much greater depth than the UV and Vis— a definite advantage in tissue imaging and therapeutic applications. In photobiology, the detection of singlet oxygen and development of efficient photo sensitizers for PDT has been the dominant application for years. The continuing introduction of new NIR emitters, coupled with better detection and lower cost systems, continues to fuel the growth of NIR luminescence applications. HORIBA offers an extensive line of NIR photoluminescence systems with a broad range of options and accessories. The detectors offered include both NIR PMTs and solid state photodiodes that span the range up to 5,500 nm. Most of these detectors can also be used with pulsed light sources, thus providing lifetime capabilities in NIR.

Applications of NIR Fluorescence


- Materials Science
 - Nanomaterials
 - Glasses and ceramics
 - LEDs and lasing media
 - Semiconductors
 - Upconverting nanoparticles for tissue imaging
- Optical fiber communication
 - Optical amplifiers (e.g. chelated Er^{3+} , 1540 nm)
- Photobiology and photomedicine
 - Singlet oxygen detection
 - R&D of singlet oxygen photosensitizers
 - Photodynamic Therapy (PDT)
- Environmental
 - Photo-oxidation of environmental pollutants
- Geology
 - NIR luminescence of minerals
- Forensic science
 - Identifying forged documents




Applications and Examples: NIR Fluorescence Spectroscopy

Photoluminescence of Lanthanides

Many applications using photoluminescence measurements involve rare earth ions (lanthanides), such as Nd^{3+} , Er^{3+} , Tm^{3+} , Ho^{3+} and Pr^{3+} , which often emit in the NIR. Often these ions are used with ligand photosensitizers which improve their light absorption properties, as lanthanide ions themselves are very weak absorbers. They are used as dopants in lasing media and glasses, and are made into nanoparticles of varying sizes and shapes in order to control their optical properties. The photoluminescence lifetime (from microseconds to milliseconds) is the key parameter in assessing the optical efficiency of devices involving lanthanides, as well as in quality control during their manufacturing.



PL emission and excitation spectra of Ho^{3+} doped glass measured with the PTI QuantaMaster 8000, using the TE-cooled InGaAs detector and lock-in amplifier.




PL decays of Ho^{3+} doped glass measured with the PTI QuantaMaster 8000 system operating in the lifetime mode. Note that the decays are very different for different transitions. The decay at 1200 nm also shows a rise time of 2.4 μs .


Single Walled Carbon Nanotubes (SWCNT)

Single-walled carbon nanotubes is a hot topic in photonics/materials science in the last few years. There are numerous existing and potential applications where SWCNTs are used, such as microdiodes and microtransistors, computing and switching devices, screen displays, gas sensors, biological sensors (NA hybridization), bio-imaging, drug delivery and many others. Mechanically, they are 3 to 10 times stronger than steel and exhibit high thermal and electrical conductivity.

SWCNTs are made of a sheet of graphene rolled along a certain angle (chiral angle) into a tube of diameter r . These structural parameters can be determined by photoluminescence measurements, usually in the NIR range. By collecting a 3D excitation-emission matrix and determining the excitation and emission wavelengths of the 3D PL peaks, the structural parameters, the chiral angle and r , can be calculated. FelixGX provides the Nanotube Calculator which makes this task easy.




Photoluminescence spectra of two different SWCNTs measured with the PTI QuantaMaster 8000 NIR photoluminescence spectrometer. In order to determine the chiral angle α and the nanotube radius r , a 3-D ExEm matrix needs to be acquired, and the result submitted to the Nanotube Calculator.



Optical Fiber Communications


The use and demand for optical fiber communication has grown rapidly, and applications are numerous, ranging from global networks to desktop computers. There have been three spectral 'windows' used for optical transmission: 850 nm, 1310 nm and 1550 nm, with the third window now becoming a globally accepted transmission band. There is a need to insert some light amplifiers along the fiber line. One idea of amplifying the signal is based on a chelated Erbium ion. Erbium belongs to the family of lanthanides and has an emission band in the NIR at about 1550 nm, so it matches perfectly the 3rd optical transmission window. The chelating molecules are excited in the UV or Vis by inexpensive LEDs and transfer the excitation energy by FRET to the Erbium center, thus promoting Erbium to its excited state. Since the energy difference between the excited and ground state of Erbium equals the energy of photons (1550 nm) that are propagated along the fiber, these incoming photons will stimulate the emission from Erbium, enhancing the overall signal.




Emission spectrum of chelated Erbium (solid sample) measured with the NIR-PMT and solid sample holder accessory. The system features a thermoelectrically cooled, extended wavelength range NIR PMT operating in the photon counting mode.

Singlet Oxygen Detection


Singlet oxygen generation and detection are growing fields with applications in such areas as cancer treatment, photosensitized oxidations, and biomolecular degradation. The first excited state of an oxygen molecule is a singlet state, which can readily react with other singlet molecules. Radiative decay to the triplet ground state is a spin forbidden transition resulting in a long lived excited state. Excited singlet oxygen emits phosphorescence in the NIR at 1270 nm.



Spectra of singlet oxygen generated by Rose Bengal in methanol



Lifetime of singlet oxygen using a NIR-PMT with the Xe flash lamp. This is a substantial undertaking, considering the singlet oxygen phosphorescence quantum yields are of the order of 10^{-6}



Luminescence decay of chelated Erbium (solid sample) measured with the NIR-PMT system operating in the time-resolved 'gated' mode. The decay can be described by a broad tri-modal lifetime distribution, as shown by the MEM distribution analysis—a powerful analysis package from HORIBA.

Accessories

Four Position Peltier K-157- C Temperature Control

Fully automated Peltier-based temperature control (-25°C to 105°C) for up to 4 samples measured at pre-set temperatures, or with temperature ramping and simultaneous measurement at all 4 sample positions. Magnetic stirring included.



K-Uni-HDLR Universal Holder for Solids, Powders and Cuvettes

Universal Sample Holder Base, capable of both linear and rotational travel, was designed for the measurements of solid compounds, microscope slides, or films. The solid sample holder head, also suitable for slides and films, mounts onto a base and can be removed easily to substitute a powder sample holder head or a cuvette holder which enables front face measurements.

Cold Finger Dewar K-158

The cold finger dewar accessory is designed to be used with liquid nitrogen as coolant (77 K). Includes: quartz cold finger dewar that accepts 5 mm tubes, dewar holder for the sample turret or single cuvette holder, foam lid for the dewar and extension collar with altered sample chamber lid, and a sample compartment. The dewar features a suprasil quartz cold finger that passes light down to about 200 nm. Samples are placed in NMR and EPR tubes, and the liquid nitrogen placed in the dewar will typically last several hours.



Titration K-165-B

Titration is performed to measure a number of biochemical and physical parameters, including binding constants, stoichiometry and kinetics. HORIBA offers fully automated titration solutions that are integrated into the software. Parameters such as mixing, volume, speed, and calibration are dictated in the software and can be adapted to your needs.

Stopped Flow Accessory K-161-B


The stopped flow accessory is used to rapidly mix small volumes of two (or more) different chemicals in a cuvette, quickly stop the flow of chemicals to the cuvette, and monitor the resulting chemical reaction via optical means. In some instances, the chemical reaction will result in luminescence, and this optical signal can be monitored using a fluorometer. In other instances, the chemical reaction only produces a change in the optical absorption properties and must be monitored using an absorption technique. The primary experimental interest is in the speed of the chemical reaction following the mixing in the cuvette, in addition to the spectral properties of the resulting absorption and/or luminescence.

Cryostat K-CRYO-3/K- CRYO-2

The Cryostat can be used between the temperatures of 77 K up to 500 K. FelixGX can control the cryostat remotely to allow computer monitoring of steady state and lifetime measurements.

Muscle Strip Accessory K-162-B

The muscle strip is inserted into a standard 1 cm cuvette, combining the lower muscle hook with unique perfusion tubes, a tension transducer with upper muscle hook, and an interface electronic control unit. The accessory can be used with any cuvette-based fluorescence system having a standard single cuvette holder complete with tension transducer and transducer mounting bracket with micrometer position adjustment.



Single Cuvette Peltier K-155-C Temperature Control

The Peltier-based temperature control with magnetic stirring provides unmatched temperature stability and full software control from -40°C to 105°C , including temperature ramping experiments.



Remote Probe Sensing Accessory K-163

The remote probe sensing accessory allows in vitro or in vivo measurement by means of a quartz bifurcated fiber bundle or Liquid Light Guide. One fiber leg is attached to the second exit port of the excitation monochromator to provide excitation light to the sample. The second leg is attached to an open entrance port of the emission monochromator to detect the fluorescence signal emitted from the sample.

Polarizers

PTI offers a wide variety of polarizers, ranging from manual sheet polarizers to automated large aperture Glan Thompson polarizers. All configurations allow for automated software control, automatic G-factor determination, and real-time acquisition of HH, VH, VV, and HV analysis. Measure steady state anisotropy in a single emission configuration or dynamic anisotropy utilizing our dual emission configuration.

Total Internal Reflection Fluorescence (TIRF) Flow System Accessory K-TIRF

Incorporate TIRF flow cell into a PTI system by replacing the standard sample holder. Includes two reusable UV silica TIRF prisms, twenty reusable silica TIRF slides, two plastic fluidic blocks, a set of ten elastic gaskets, and a Teflon cassette holder for cleaning and chemical modification of ten TIRF slides.


Absorption Accessory ABS-ACC

The ABS-ACC absorption accessory fits directly into a cuvette sample holder and enables the user to measure the absorption spectrum, or check the optical density of the sample without reconfiguring the PTI QuantaMaster 8000 fluorometer.

Integrating Sphere K SPHERE-B

High Performance Integrating Sphere: Designed for enhanced measurement of quantum yields of solids, films, and powders.

We use a 6-inch diameter sphere and attach it directly to the sample chamber on the port opposite the excitation channel. This design minimizes the effect of the excitation, emission, and sampling ports on the accuracy of the measurement.




Integrating Sphere K SPHERE-Petite

Petite Integrating Sphere: Easy to use 3.2-inch Integrating Sphere that can be installed in seconds, replacing a standard cuvette holder - no optics to change or light guide coupling required. Easily removable top enables changing samples in seconds. It can accommodate regular cuvettes, slides and powders, with all 3 holders included.

HORIBA FLUORESCENCE

Scientific



www.horiba.com/fluorescence
info.sci@horiba.com

- USA:** HORIBA Instruments Inc., 3880 Park Avenue, Edison, NJ 08820-3012 - Toll-free: +1-866-562-4698 - Tel: +1 732 494 8660 - Fax: +1 732 549 5125 - Email: info.sci@horiba.com
- France:** HORIBA Jobin Yvon S.A.S., 16-18 rue du Canal, 91165 Longjumeau cedex - Tel: +33 (0)1 69 74 72 00 - Fax: +33 (0)1 69 09 07 21 - Email: info-sci.fr@horiba.com
- Japan:** HORIBA Ltd., Tokyo Branch Office, 2-6, KandaAwaji-cho, Chiyoda-ku, Tokyo 101-0063, Japan - Tel: +81-(0)3 6206 4721 - Fax: +81 (0)3 6206 4730 - Email: info-sci.jp@horiba.com
- Germany:** HORIBA Jobin Yvon GmbH, Hauptstrasse 1, 82008 Unterhaching - Tel: +49 (0)89 4623 17-0 - Fax: +49 (0)89 4623 17-99 - Email: info-sci.de@horiba.com
- Italy:** HORIBA Jobin Yvon Srl, Via Cesare Pavese 21, 20090 Opera (Milano) - Tel: +39 2 5760 3050 - Fax: +39 2 5760 0876 - Email: info-sci.it@horiba.com
- UK:** HORIBA UK Ltd., 2 Dalston Gardens, Stanmore, Middlesex HA7 1BQ - Tel: +44 (0)20 8204 8142 - Fax: +44 (0)20 8204 6142 - Email: info-sci.uk@horiba.com
- China:** HORIBA (China) Trading Co. Ltd., Unit D 1F, Bldg A, Srynnex International Park, No. 1068 West Tianshan Road, Shanghai 200335 - Tel: +86 (0)21 6289 6060 - Fax: +86 (0)21 6289 5553 - Email: info-sci.cn@horiba.com
- Brazil:** HORIBA Instruments Brasil Ltda., Rua Presbítero Plínio Alves de Souza, 645, Loteamento Polo Multivias, Bairro Medeiros, Jundiaí / SP, CEP 13.212-181 - Tel: +55 (0)11 2923 5400 - Fax: +55 (0)11 2923 5490 - Email: info.ciencia.br@horiba.com
- Other:** Tel: +1 732 494 8660 - Email: info.sci@horiba.com

HORIBA Scientific has a policy of continuous product development, and reserves the right to amend part numbers, descriptions and specifications without prior notice.