

TGS-Series

Toroidal Grating Spectrograph

Coming from Synchrotron technologies

Compact

Stable

True Flat Field

No Motor

Robust

Fixed position

High Throughput

Grazing Incidence




A spectrograph to explore the 9.5–170 nm spectral range

Variable Line Spacing Gratings

TGS use aberration-corrected VLS toroidal gratings which disperse and refocus the light from the entrance slit onto the exit focal plane of the spectrograph.

The groove spacing of these gratings is computer-optimized to produce high quality images with a minimum of astigmatism and coma over a large spectral range and even at high numerical aperture.

The VLS grating grooves are no longer straight and parallel, but instead correspond to confocal hyperboloids or ellipsoids. Optimizing the position, angles and arm lengths of the two recording beams provides the optical designer with the degrees of freedom necessary to minimize aberrations.



Two gratings in slider

The TGS (Toroidal Grating Spectrograph) series is especially designed for analyzing the whole VUV wavelength range from 9.5 to 170 nm with high spectral resolution thanks to its bi-dimensional flat field.



TGS300 in UHV version, with 2 gratings on slider

A full spectrum range is simultaneously acquired in a single acquisition when adding an array detector. With a standard DN100CF exit flange mounted on XY table, a 25 mm array detector will be easily optimized in wavelength centering and focusing over the 40 mm flat spectrum.

Applications

- **High Harmonic Generation**
- **EUV and UV Plasma analysis**
- **Light source characterization**
- **VUV Laser Analysis**

Features

- Single Toroidal Grating design
- Choice of three gratings
- Low astigmatism level
- Fixed grating position
- CCD flange adjustable in exit focal plane
- Choice of master or replica gratings
- Kinematic grating mount or two gratings slider
- Compact design
- High Vacuum compatible


Benefits

- Optimized for throughput
- Optimized spectral range and resolution
- 40 mm Flat field spectrograph
- Robust
- In vacuum focus and central wavelength adjustment
- Optimize damage threshold or cost
- Interchangeable gratings
- Less room consuming and easy to install
- HV (10^{-6} mbar) – Optional UHV operation (10^{-9} mbar)



A fixed system with flat field focal plane

A grazing VLS designed to optimize resolution and efficiency

Based on HORIBA Scientific technology, the TGS series is designed around a single toroidal aberration corrected grating working in grazing incidence to increase the efficiency in the EUV range. Instead of having a fixed spherical grating and the exit port rotating on the Rowland circle as most of the other common design, the TGS has the particularity to have a flat field focal plane almost perpendicular to the exit axis of the instrument. Therefore, the CCD detector is fixed in position and acquires simultaneously a full spectrum with an unmatched throughput and resolution over the whole spectral range, making the TGS a very compact and stable solution.



TGS300 optical layout



High harmonic generation from multiple orbitals in N_2 (B. K. McFarland, J. P. Farrell, P. H. Bucksbaum and M. Gühr, *Science*, **322**, 1232 (2008))

Pinch discharge operated in different working gases recorded with a custom TGS (HEXOS 1)

The calibration of the instrument with harmonics and plasma emission is described in: Calibration of a High Harmonic Spectrometer by Laser Induced Plasma Emission, J. P. Farrell, B. K. McFarland, P. H. Bucksbaum and M. Gühr, *Optics Express*, **17**, 15134-15144 (2009)

Gratings available

Model	Gratings		Spectral range		Exit Dispersion	Resolution (FWHM in nm)
	Part Number	Density (gr/mm)	nm	eV		
TGS300	541 00 220	2105	9.5 - 32	39 - 130	0.5 nm/mm at 9.5 nm, 0.6 nm/mm at 32 nm	2105 gr/mm 0.036 nm
	541 00 200	450	10 - 110	11 - 124	2.3 nm/mm at 10 nm, 2.7 nm/mm at 110 nm	450 gr/mm 0.075 nm
	541 00 210*	290	15.5 - 170	7 - 80	3.5 nm/mm at 15.5 nm, 4.2 nm/mm at 170 nm	290 gr/mm 0.11nm

* Gratings are available in master or replica version


Options

- Ultra High Vacuum (UHV, 10^{-9} mbar) version
- Slider for grating change under vacuum
- Removable entrance arm
- Laser kit for easy alignment




Accessories

- EUV/UV light source
- CCD detectors




	TGS300
Optical design	Fixed toroidal VLS Grating (single optic)
Focal length	300 mm
Aperture	f/30
Optic coating	Pt (Au in option)
Deviation angle	137°
Vacuum	10 ⁻⁶ mbar (HV version) or 10 ⁻⁹ mbar (UHV version)
Pumping flange	DN63LF (DN63CF)
Entrance port	Micrometric slit (10 μm to 2 mm)
Exit port	Adjustable CCD port
Entrance flange	DN40 KF
Exit flange	DN100CF
Software	HORIBA Scientific software
PC Interface	RS232 - USB2



Theoretical Absolute Efficiency



TGS300 in HV version, two gratings on slider



TGS300 in UHV version, two gratings on slider



www.horiba.com/scientific
info.sci@horiba.com

France: HORIBA Jobin Yvon S.A.S., 16-18 rue du Canal, 91165 Longjumeau cedex - Tel: +33 (0)1 69 74 72 00 - Fax: +33 (0)1 69 09 07 21 - Email: info-sci.fr@horiba.com
USA: HORIBA Instruments Inc., 3880 Park Avenue, Edison, NJ 08820-3012 - Toll-free: +1-866-562-4698 - Tel: +1 732 494 8660 - Fax: +1 732 549 5125 - Email: info-sci.us@horiba.com
Japan: HORIBA Ltd., Tokyo Branch Office, 2-6, KandaAwaji-cho, Chiyoda-ku, Tokyo 101-0063, Japan - Tel: +81-(0)3 6206 4721 - Fax: +81 (0)3 6206 4730 - Email: info-sci.jp@horiba.com
Germany: HORIBA Jobin Yvon GmbH, Hauptstrasse 1, 82008 Unterhaching - Tel: +49 (0)89 4623 17-0 - Fax: +49 (0)89 4623 17-99 - Email: info-sci.de@horiba.com
Italy: HORIBA Jobin Yvon Srl., Via Cesare Pavese 21, 20090 Opera (Milano) - Tel: +39 2 5760 3050 - Fax: +39 2 5760 0876 - Email: info-sci.it@horiba.com
UK: HORIBA UK Ltd., 2 Dalston Gardens, Stanmore, Middlesex HA7 1BQ - Tel: +44 (0)20 8204 8142 - Fax: +44 (0)20 8204 6142 - Email: info-sci.uk@horiba.com
China: HORIBA (China) Trading Co. Ltd., Unit D 1F, Bldg A, Srynnex International Park, No. 1068 West Tianshan Road, Shanghai 200335 - Tel: +86 (0)21 6289 6060 - Fax: +86 (0)21 6289 5553 - Email: info-sci.cn@horiba.com
Brazil: HORIBA Instruments Brasil Ltda., Rua Presbítero Plínio Alves de Souza, 645, Loteamento Polo Multívias, Bairro Medeiros, Jundiaí / SP, CEP 13.212-181 - Tel: +55 (0)11 2923 5400 - Fax: +55 (0)11 2923 5490 - Email: infocientifica.br@horiba.com
Other: Tel: +33 (0)1 69 74 72 00 - Email: info.sci@horiba.com

HORIBA

Scientific